ManuStabil Handgelenkorthese

ManuStabil Wrist Brace

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ManuStabil Wrist Brace

- With reinforcement splints
- Breathable material

Many thanks for placing your trust in this medical device. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

INTENDED PURPOSE

This medical device is a wrist brace for immobilisation in at least two directions of movement.

INDICATIONS

Pre-operative, post-operative, post-traumatic, degenerative, carpal tunnel syndrome (CTS), tenosynovitis, peripheral nerve palsy.

CONTRAINDICATIONS

Peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders with the region of the body treated, skin diseases in the part of the body treated.

APPLICATION RISKS/IMPORTANT NOTES

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- do not close it too tightly during longer periods of rest, e.g. when sleeping
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only

- whilst wearing the brace, please neither use any creams nor ointments in the vicinity of the attached aid as they can destroy the material

FITTING INSTRUCTIONS

The brace can be used on both sides. Here, ensure that the thumb is either guided through the left or the right thumb opening, dependent on the side. Before putting on, bend the aluminium splints to the desired shape. Guide the hand into the brace and firstly set the metacarpus circumference on the back of the hand using the two Velcro fasteners. Then, close the Velcro fastener on the lower arm. Finally, attach the circumferential strap to the wrist.

Tip: The metacarpus strap can be left in the guide loops to enable easier removal and fitting.

REMOVAL

To remove the brace, please open the Velcro fasteners. Grasp the splint in the palm area and pull the brace off over the fingers. Close the Velcro fasteners in order to avoid early wear to the Velcro surface.

MATERIAL COMPOSITION

Polyamide (PA), polyester (PES), polyurethane (PUR)

The sewn in textile label provides the precise material composition.

CLEANING INFORMATION

 Delicate wash 30 °C  Do not bleach  Do not dry in a tumble dryer

 Do not iron  Do not clean chemically

Remove the aluminium splint before washing. Close the Velcro fasteners in order to avoid damaging other items of clothing. Do not use fabric conditioner. Stretch back into shape and dry in the fresh air.

GUARANTEE

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

USEFUL LIFE/LIFETIME OF THE PRODUCT

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

DUTY OF NOTIFICATION

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find the manufacturer's contact data in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/md-eu-contact.

DISPOSAL

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

DECLARATION OF CONFORMITY

The manufacturer confirms that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL.

Status: 04.2023

 Medical device |  Single patient – multiple use

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VPS Medical Ltd.
BORT GmbH | Am Schweizerbach 1 | D-71384 Weinstadt