COMFORT Handgelenkbandage

COMFORT Wrist Support

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COMFORT Wrist Support

- SOFTflex Plus knit fabric
- Cooling luxury fibre

Many thanks for placing your trust in this medical device. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

INTENDED PURPOSE

This medical device is a wrist compression support for the compression of soft tissue with a local pressure pad and a reinforcement splint.

INDICATIONS

Chronic, post-traumatic or post-operative irritations in the wrist area, e.g. distortion, carpal arthrosis, tenosynovitis.

CONTRAINDICATIONS

Peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders with the region of the body treated, skin diseases in the part of the body treated, latex allergy.

APPLICATION RISKS/IMPORTANT NOTES

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- remove it during longer periods of rest, e.g. when sleeping
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used

- no re-use. This medical aid is intended for treating one patient only
- whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material
- In order to avoid damages to the knitted support fabric, only attach the circular stabilizing strap after successfully determining the size and trial fitting. Goods damaged by the Velcro fastening cannot be returned.

FITTING INSTRUCTIONS

Place the support on the metacarpophalangeal joint affected and guide the thumb through the intended thumb opening at the same time. The plastic splint runs along the surface of the palm. The circular stabilizing fleece strap included can be fitted individually dosed. For this, attach the circular stabilizing fleece strap to the velour on the support intended for this purpose using the narrow Velcro strip. The circular stabilizing fleece strap can then be twisted around the wrist circumferentially. If required, the plastic splint material can be individually adapted by an orthopaedic technician at a temperature of approx. 180°C. Let the reshaped splint cool down and then re-insert it into the support splint bag.

REMOVAL

To remove the support, please loosen the circular stabilizing fleece strap. Grasp the splint in the palm area and pull the support off over the fingers. Close the Velcro fastener in order to avoid early wear to the Velcro surface.

MATERIAL COMPOSITION

Polyamide (PA), elastodiene/latex (ED), elastane (EL)

The sewn in textile label provides the precise material composition.

The product contains latex and can trigger allergic reactions.

CLEANING INFORMATION

Delicate wash 30°C Do not bleach Do not dry in a tumble dryer

Do not iron Do not clean chemically

Remove the plastic splint before washing. Close the Velcro fastener in order to avoid damaging other items of clothing. Do not use fabric conditioner. Stretch back into shape and dry in the fresh air.

GUARANTEE

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

USEFUL LIFE/LIFETIME OF THE PRODUCT

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

DUTY OF NOTIFICATION

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find the manufacturer's contact data in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/md-eu-contact.

DISPOSAL

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

DECLARATION OF CONFORMITY

The manufacturer confirms that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL.

Status: 09.2023

 Medical device |  Single patient – multiple use

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VPS Medical Ltd.
BORT GmbH | Am Schweizerbach 1 | D-71384 Weinstadt